

MAN IS A WONDER
CONSISTING OF
80 TRILLION CELLS.
MAN IS A
BIOMAGNET.

A person's back is shown from the waist up, wearing a white t-shirt. A large, circular, blue device with a black border is placed on the lower back. The device has a subtle grid pattern and the text 'epadR' in white. A handheld probe, which is black and silver, is connected to the device by a coiled black cord. The probe is held against the person's skin. The background is a solid dark blue.

epad^R



MOBILE ENERGY REFUELLING STATION.



RECHARGING YOUR BATTERIES IS NOW POSSIBLE WHEREVER YOU ARE. THE NEW MOBILE ePADS OFFER YOU A SIMPLE WAY TO BOOST THE IMMUNE SYSTEM, ACCELERATE REGENERATION AND RECHARGE OUR BATTERIES.



CARPE DIEM

The use of magnetic field therapy has never been easier. The ePAD developers have worked extensively to make this form of cell renewal easily mobile. You can now recharge your batteries while travelling in a car or on a plane- a perfect opportunity to use the extra time to improve your energy balance – in a simple

and flexible way. You can even refuel while doing household chores, watching TV or relaxing at home.

Three different programmes allow you to adjust the device to your personal needs. Small aches and pains can also easily be integrated into regeneration.



PROGRAMME SELECTION

PROGRAMME No. 1 (8 Hertz)

PRESS THE BLUE BUTTON ONCE: The two waves start to flash red every 4 seconds. This provides the same frequency we have in our normal waking state and provides support for the body.

PROGRAMME No. 2 (5 Hertz)

PRESS THE BLUE BUTTON TWICE: The two waves start to flash blue every 4 seconds. This provides the same frequency the body has during sleep and provides the body with relaxation.

PROGRAMME No. 3 (3 Hertz):

PRESSING THE BLUE BUTTON THREE TIMES: The two waves start to flash green every 4 seconds. This provides the same frequency the body has in deep sleep, when the body completely relaxes and tensions are released.

PRESSING THE BLUE BUTTON FOUR TIMES: The device is switched off. The device is in operation only when both waves flash at the same time every 4 seconds. Each programme runs for 25 minutes and then the device turns off automatically and the programme is complete.



USE THE STRENGTH OF THE EARTH IN SPORTS.



AFTER EXERCISE, IT IS ESSENTIAL TO REGENERATE QUICKLY AND REGAIN ENERGY.



CARPE NATURA

The primary aim of sport and exercise is to draw on the body's numerous energy sources. Those sources can become fatigued, resulting in pain and sometimes injury. The ePAD Sport can promote optimal energy flow through better blood circulation, the release of certain blockages and relieve muscle cramping or soreness

– all while on the go – during car and train rides, flights, at work, in front of the TV, reading or relaxing.

AREAS OF APPLICATION:

After exercise-oxygen supply can reach the muscles for faster recovery time. The gentle stimulation also supports the healing process of acute injuries.

LEGS – KNEES – ARMS – BACK – NECK

PROGRAMME SELECTION

PROGRAMME No. 1 (1500 Hertz)

PRESSING THE GREEN BUTTON ONCE:

The two waves start to flash red every 4 seconds. This frequency supports acute injuries and has an anti-inflammatory effect.

PROGRAMME No. 2 (8 Hertz)

PRESSING THE GREEN BUTTON ONCE:

The two waves start to flash blue every 4 seconds. This is similar to body in its active waking state, providing energetic support for the body.

PROGRAMME No. 3 (3 Hertz):

PRESSING THE BLUE BUTTON THREE

TIMES: The two waves start to flash green every 4 seconds. This provides the same frequency the body has in deep sleep, when the body completely relaxes and tensions are released.

PRESSING THE BLUE BUTTON FOUR

TIMES: The device is switched off. The device is in operation only when both waves flash at the same time every 4 seconds. Each programme runs for 25 minutes and then the device turns off automatically and the programme is complete.





THE FOCUS IS ON FINDING
BALANCE OF BODY AND MIND.
THIS IS THE ONLY WAY FOR
THE BODY TO REGENERATE.